Biosan SIA, Ratsupites iela 7 k-2, Riga, Latvia, LV-1067 E-mail: marketing@biosan.lv, www.biosan.lv Phone:+37167426137

SC-18, block

DESCRIPTION

 20×0.5 ml + 12×1.5 ml microtubes

SPECIFICATIONS

Amount of 0.5 ml tubes	20
Amount of 1.5 ml tubes	12
Max. RPM	1400
Max. temperature	100°C



SC-18

CAT. NUMBER

BS-010120-AK







TS-100 Without thermoblock Thermo-Shaker for microtubes and PCR plates

Thermo–Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

read more