Biosan SIA, Ratsupites iela 7 k-2, Riga, Latvia, LV-1067 E-mail: marketing@biosan.lv, www.biosan.lv Phone:+37167426137



VP-4, block

DESCRIPTION

 4×50 ml conical tubes

SPECIFICATIONS

			Statement of the owner owner owner owner owner
Max. RPM		1000	
Max. temperature	8	30°C	



VP-4

CAT. NUMBER

BS-010175-GK





TS-100 Without thermoblock Thermo-Shaker for microtubes and PCR plates

Thermo–Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

read more