

Features:

- Highly resistant towards a broad range of chemicals
- Easy selection and dispensing through a vision panel
- Reliable and maintenance-free due to its robust design
- Flexible: 12 tip sizes, 10 volume adjustments and 120 program steps
- Reliable ergonomics for your lab routine without sacrificing ultimate precision
- · Handy, having a perfect ergonomic design. One-handed operation of volume selection, loading and dispensing of liquids
- Convenient: the easy-moving dosing button and the lightweight instrument makes long dosing series more comfortable
- The CAPP RIO Mechanical Repeater is highly compatible with both CAPP Harmony syringes and also with tips of many common brands
- The table below shows all dosing possibilities, enabling you to choose the most suitable tip for your application



CAPP % Dispensing Volumes

Program No.	1	2	3	4	5	6	7	8	9	10
0.05mL	1μL	1.5µL	2μL	2.5µL	3μL	3.5µL	4µL	4.5μL	5μL	5.5μL
0.1mL	2μL	3μL	4μL	5μL	6μL	7μL	8μL	9μL	10μL	11μL
0.2mL	4μL	6μL	8μL	10μL	12μL	14μL	16μL	18μL	20μL	22μL
0.5mL	10μL	15μL	20μL	25μL	30μL	35μL	40μL	45μL	50μL	55μL
1mL	20μL	30μL	40μL	50μL	60μL	70μL	80μL	90μL	100μL	110µL
1.25mL	25μL	37.5μL	50μL	62.5μL	75μL	87.5μL	100μL	112.5μL	125μL	137.5μL
2.5mL	50μL	75μL	100μL	125μL	150μL	175μL	200μL	225μL	250μL	275μL
5mL	100μL	150μL	200μL	250μL	300μL	350μL	400μL	450μL	500μL	550μL
10mL	200μL	300μL	400μL	500μL	600μL	700μL	800μL	900μL	1000μL	1100μL
12.5mL	250μL	375μL	500μL	625μL	750μL	875μL	1000μL	1125μL	1250μL	1375μL
25mL	500μL	750μL	1000μL	1250μL	1500μL	1750μL	2000μL	2250μL	2500μL	2750μL
50mL	1000μL	1500μL	2000μL	2500μL	3000μL	3500μL	4000μL	4500μL	5000μL	5500μL
Steps	48	31	22	18	14	12	10	9	8	7

Ordering Information

Cat. No.	Description				
PR-10	CAPP R10 Mechanical Repeater				

